



## Individual Sieve Analysis

Finding the right sieve shaker is easy: Simply send us a sample of your choice – we will conduct a sieve analysis and send you an individual sieving report and recommend an instrument suitable for your application.

Please complete the form completely and **email it in advance to [deoliveira@fritsch-france.fr](mailto:deoliveira@fritsch-france.fr) and send us the material together with the print out of the completed form.**

If you would like to send an additional sample (max. 2 samples) which differs in regards to consistency, desired sample quantity or final fineness, please complete a second form for this second sample.

**The fields marked with an asterisk\* are required fields and have to be completed!**

### Your information about the material

Name of the material\*:

Chemical formula:

**Hazard material\*:**  yes<sup>1</sup>  no

(\*Please enclose safety data sheet!)

explosive  toxic  caustic  oxidising  environmental hazard

easily flammable  harmful to health from:

May not be put in contact with

### Material properties

hygroscopic humidity:  %

The material may be dried / heated up to:  °C

Soluble in:

Other:

### Task

Which quantity should be sieved per charge \*:  g

\* depends on sample and utilized sieves

<b>Dry sieving</b>	Vibratory Sieve Shaker ANALYSETTE 3 Heavy Duty Analytical Sieve Shaker ANALYSETTE 18	for sieves < 63 mm: up to 2 kg*, for sieves < 100 µm: up to 100 g* up to 15 kg*
<b>Wet sieving</b>	Vibratory Sieve Shaker ANALYSETTE 3 Heavy Duty Analytical Sieve Shaker ANALYSETTE 18	20 - 100 g* up to 1 kg*
<b>Micro precision sieving</b>	Vibratory Sieve Shaker ANALYSETTE 3 PRO	0.05 - 0.5 g* with max. 4 micro precision sieves

What type of sieving do you request?

Dry sieving  Wet sieving  Micro precision sieving

Which sieving aids may be used with dry sieving?

none  agate balls 5 / 10 mm  rubber balls 20 mm

Vulcollan cubes  dispersing agent

May wetting agents with wet sieving in water be used?

Yes, we recommend:   no

What kind of liquid do you recommend for micro-precision sieving?

water  other:

### Which sieve shaker should be utilized?

- Please select the suitable instrument for our task!
- Vibratory-Sieve Shaker ANALYSETTE 3 PRO
- Vibratory-Sieve Shaker ANALYSETTE 3 SPARTAN
- Heavy Duty Analytical Sieve Shaker ANALYSETTE 18



**Sieves with the following mesh widths should be used for the sieve analysis?\***

<input type="checkbox"/> mm	<input type="checkbox"/> µm	<input type="checkbox"/> mesh
a) <input type="text"/>	b) <input type="text"/>	c) <input type="text"/>
d) <input type="text"/>	e) <input type="text"/>	f) <input type="text"/>
g) <input type="text"/>	h) <input type="text"/>	i) <input type="text"/>

How did you conduct the particle size analysis in the past?

**Which results did you obtain?**

<input type="checkbox"/> mm	<input type="checkbox"/> µm	<input type="checkbox"/> mesh	
Aperture	Cumulative weight undersize	Aperture	Cumulative weight undersize
a) <input type="text"/> = <input type="text"/> %		b) <input type="text"/> = <input type="text"/> %	
c) <input type="text"/> = <input type="text"/> %		d) <input type="text"/> = <input type="text"/> %	
e) <input type="text"/> = <input type="text"/> %		f) <input type="text"/> = <input type="text"/> %	
g) <input type="text"/> = <input type="text"/> %		h) <input type="text"/> = <input type="text"/> %	
i) <input type="text"/> = <input type="text"/> %		j) <input type="text"/> = <input type="text"/> %	

Remarks

Would you like to receive an offer?  yes  no

Should not needed material be returned?  yes  no

**Your personal information**

Salutation*:	<input type="text"/>	Title:	<input type="text"/>
Last Name*:	<input type="text"/>	First name:	<input type="text"/>
Company*:	<input type="text"/> Please supply end customer info	Department:	<input type="text"/>
Street*:	<input type="text"/>	House No.:	<input type="text"/>
Postcode*:	<input type="text"/>	City*:	<input type="text"/>
Country*:	<input type="text"/>	E-Mail*:	<input type="text"/>
Phone*:	<input type="text"/>	Fax:	<input type="text"/>

\*Yes, I read the Privacy Policy and consent to that data supplied by me, is electronically processed and saved. My data is used exclusively for this purpose.

I consent to, that my aforementioned data is saved and used for the mailing of further information about your products, services and events. There will be no disclosure to third parties. I can revoke this consent at any time via e-mail to [info@fritsch.de](mailto:info@fritsch.de), per letter or via clicking the unsubscribe link contained in the e-mails.



---

**Please send the completed form in advance to [deoliveira@fritsch-france.fr](mailto:deoliveira@fritsch-france.fr) and send the sample material together with the print out to:**

FRITSCH GmbH • Broyage et Granulométrie  
Monsieur Walter de Oliveira • chef produit France  
29 Chemin de la Butte  
91160 Champlan • France

Portable +33 6 60 23 89 94  
Tél./Fax +33 1 69 09 72 27  
[deoliveira@fritsch-france.fr](mailto:deoliveira@fritsch-france.fr)  
[www.fritsch-france.fr](http://www.fritsch-france.fr)

Or to our headquarters in Idar-Oberstein • Germany, the complete address can be found >> [here](#)